

TECHNICAL DATA SHEET

IAS Canine Kit 450mL

DK4BWF1

TDS-DK4BWF1-v1.0



HS CODE:

Transfusion and infusion apparatus

USAGE:

Veterinary use only

INTENDED USE:

Blood bag intended for the collection of 450ml +/- 10% whole blood. The multi bag system allows closed system preparation and subsequent storage of blood components including packed red cells, plasma products (including fresh frozen plasma, cryoprecipitate and cryosupernatant) and platelets.

TECHNICAL SPECIFICATIONS:

CPD anticoagulant solution and SAGM conservative solution to ensure the quality of blood components over time.

Long collection line allows tailored gravity flow rates by modifying table height allowing rapid draw times for large donors and flexible location of suction units away from the donor. Less likelihood of overdraw due to drag caused by short collection lines.

16G siliconized sharp needle to ensure smooth venipuncture.

No snap valves in collection lines reduces the chance of clots or sluggish flow within the line.

Removal of vacutainer sampling ports rarely used by veterinary blood banks removes the chance of extra plastics getting in the way of donation.

PACKAGING:

Blood kits are provided singly-wrapped and sterilised by steam.

30 Canine Kits are contained in one box.

SHELF LIFE

2 years

TECHNICAL DATA SHEET

IAS Canine Kit 450mL

DK4BWF1

TDS-DK4BWF1-v1.0



IAS
Innovative Animal Supplies

DRAWING:

