

TECHNICAL DATA SHEET

IAS Feline Kit

CK2BWF1

TDS-CK2BWF1-v1.0



HS CODE:

Transfusion and infusion apparatus

USAGE:

Veterinary use only

INTENDED USE:

First intended use: 100ml blood collection intended for the collection of up to 100ml +/- 10% whole blood to be used as fresh whole blood or stored whole blood in CPD.

Second intended use: Full multi bag system allows the collection, closed processing and storage of 50ml +/- 10% of whole blood into packed red cells in nutrient solution (SAGM) , and plasma products.

TECHNICAL SPECIFICATIONS:

Collection and Storage of Whole Blood up to 100mL with integrated anticoagulant solution (CPD). Centrifugation and separation of whole blood into packed red cells in integrated nutrient solution (SAGM) and plasma.

Sharp 20G needles to ensure a smooth venipuncture.

Reduce volume loss by length tubing optimization.

PACKAGING:

20G needles are provided singly-wrapped and sterilised by Ethylene Oxide.

Blood kits are provided singly-wrapped and sterilised by steam.

40 Feline Kits are contained in one box.

SHELF LIFE

2 years

TECHNICAL DATA SHEET



IAS
Innovative Animal Supplies

IAS Feline Kit

CK2BWF1

TDS-CK2BWF1-v1.0

DRAWING:

20G needle:



Feline kit:

